



Disposable-Flexible Reamers Kit


Single-use flexible reamers and guide pins for ACL/PCL reconstruction procedures

- On-target tunnel and socket placement with flexible reamers and guide pin
- For single- and double-bundle reconstruction

Single-Use Flexible Reamer

Reduces risk of damage to medial femoral condyle in medial portal procedures

- Tip structure designed for minimal friction and maximum reaming accuracy
- Flexible links for easy maneuvering and optimal tunnel placement
- Cannulated for use with the Armadillo kit's flexible guide pin
- Complete size range (5-12mm, in 0.5mm increments)
- Stainless steel – for quality, reliability, safety



Single-Use Flexible Guide Pin

Flexible and sturdy for the surgeon to achieve optimal ACL tunnel location and length

- Pin architecture ensures safe and accurate bone entry and exit to create the surgeon's intended tunnel/socket
- Nickel titanium alloy body for best flexibility and steerability
- Serves as a suture passer upon completion of the drilling procedure




T.A.G. Medical Products is a leader in the manufacturing of surgical instruments, including high-quality, innovative medical solutions for arthroscopic procedures

We tag Care

ARMADRILL™ Disposable Flexible Reamers Kit

A fully disposable kit, uniquely designed with flexible reamers and guide pin for accurate and safe single- and double-bundle ACL/PCL reconstruction

- Enables the surgeon to target the optimal site on the lateral condyle for successful anatomic placement
- Reduces risk of damage to medial femoral condyle in medial portal procedures
- Flexible and sturdy for the surgeon to achieve optimal ACL tunnel location and length
- Drilling in the medial portal does not require knee hyperflexion
- Perpendicular knee flexion enables a better view of the surgical site



T.A.G. Medical Products Corporation Ltd

Kibbutz Gaaton 2513000 Israel | Tel: +972-(0)4-9858400 | www.tag-med.com